

Technical note:

Fungal staining tools to study the interaction between the beneficial endophyte *Piriformospora indica* with *Arabidopsis thaliana* roots

Khabat Vahabi¹, Joy Michal Johnson¹, Corinna Drzewiecki² and Ralf Oelmüller^{1*}

¹Institute of General Botany and Plant Physiology, Friedrich-Schiller University Jena, Dornburger Str. 159, D-07743 Jena, Germany; ²JenaBios GmbH, Orlaweg 2, D-07743 Jena, Germany; * Correspondence to: b7oera@hotmail.de

Abstract

Piriformospora indica, a primitive endophytic fungus of Sebaciales, colonizes the roots of all plant species investigated so far and promotes plant growth. It also increases seed production, uptake of nutrients and resistance to biotic and abiotic stress. Several studies have demonstrated the degree of root colonization is crucial for the benefits of both symbionts. Furthermore, since the fungus grows inter- and intra-cellularly, efficient tools are required for the detection of fungal hyphae and spores in and around the roots. Here we describe different staining methods to study the histology of *P. indica*-*Arabidopsis thaliana* interaction.

Journal of Endocytobiosis and Cell Research (2011) 77-88

Category: Technical note

Keywords: *Piriformospora indica*, fungal staining tools, plant-microbe interaction

Received: 24 April 2011; Accepted: 01 July 2011

Introduction

Piriformospora indica is an endophytic fungus of Sebaciales that colonizes the roots of many plant species including forestry, agricultural, horticultural and medicinal species (Barazani et al. 2005; Glen et al. 2002; Kaldorf et al. 2005; Oelmüller et al. 2005; Peškan-Berghöfer et al. 2004; Pham et al. 2004; Sahay and Varma 1999; Selosse et al. 2002a, b; Shahollari et al. 2005, 2007; Sherameti et al. 2005, 2008a, b; Urban et al. 2003; Varma et al. 1999, 2001; Waller et al. 2005; Weiss et al. 2004), as well as the model plant *Arabidopsis thaliana* (cf. Peškan-Berghöfer et al. 2004; Oelmüller et al. 2005; Shahollari et al. 2005; Sherameti et al. 2008a; Camehl et al. 2011). For all plant species tested so far, *P. indica* showed to have the potential to promote growth and biomass production indicating that the symbiosis is beneficial for the hosts and has the potential for agricultural and biotechnological applications. We use *Arabidopsis* as host to understand the molecular basis of the beneficial interaction between the symbionts.

P. indica hyphae can be detected on the surface of Ara-

bidopsis root, in the outer cell layers, and within the root cells, using optical microscopy (Peškan-Berghöfer et al. 2004). Several studies have demonstrated that the degree of root colonization is critical for establishing benefits for the plants. For the wild-type, root colonization and association of fungal hyphae results in the promotion of growth and seed yield, as well as tolerance against biotic and abiotic stress (Varma et al. 1999; Verma et al. 1998; Oelmüller et al. 2009; Oelmüller et al. 2004; Rai et al. 2004; Sherameti et al. 2008a, b; Varma et al. 2001; Waller et al. 2005), whereas overcolonization of hyphal growth can result in the loss of the benefits, and may even lead to a shift from a mutualistic to a antagonistic interaction (Johnson and Oelmüller 2009; Camehl et al. 2010a, b). Furthermore, it has been proposed that the beneficial interaction of *P. indica* with barley roots requires cell death programs (Deshmukh et al. 2006; Schäfer et al. 2007). Fungal colonization increases with root tissue maturation. The root tip meristem of barley showed no colonization, and the elongation zone showed mainly intercellular colonization. In contrast, the differentiation zone was heavily infested by inter- and intracellular hyphae and intracellular chlamydospores. The majority of hyphae were present in dead rhizodermal and cortical cells that became completely filled with chlamydospores. In some cases, hyphae penetrated cells and built a meshwork around plasmolyzed protoplasts, suggesting that the fungus either actively kills cells or senses cells undergoing endogenous programmed cell death. Further studies with the barley BAX inhibitor-1, a gene capable of inhibiting plant cell death, demonstrated that *P. indica* requires host cell death for proliferation in differentiated barley roots (Deshmukh et al. 2006; Schäfer et al. 2007). Recently, this study was extended to *Arabidopsis* (Jacobs et al. 2011). By employing transmission electron microscopy and epifluorescence microscopy along with reporter and mutant plants, the authors observed a biotrophic colonization of *Arabidopsis* roots by a fungus. This biotrophic stage is followed by a cell death-associated colonization phase. The genetic and molecular analyses also demonstrated the efficiency of the root innate immune system to halt microbial colonization and indicated that mutualistic colonization success is intimately dependent on efficient immune suppression strategies (Jacobs et al. 2011).

All these observations clearly demonstrate that powerful staining techniques are required to visualize hyphal growth within the roots.

Until now, the growth and morphology of *P. indica* hyphae and spores was studied by light microscopy, epifluorescence, scanning-electron (Sun et al. 2010) or confocal microscopy (Stein et al. 2008; Waller et al. 2005). However,

some of the light microscopical techniques were time-consuming and harsh (Verma et al. 1998). Recently, GFP-transformed *P. indica* strains have been reported (Zuccaro et al. 2009), which might be a powerful tool in future studies.

We use *Arabidopsis thaliana* as a host for *P. indica* and employed various simple methods for the visualization of the spores and hyphae within and around the root. This includes live cell microscopy, different staining methods and epifluorescence microscopy. The objectives of this study were to develop fast and simple methods to describe the colonization under different growth and cultivation conditions, and to obtain staining protocols that allow study of the fungal growth to some extent in and around the *Arabidopsis* roots.

Results and Discussion

The interaction of *P. indica* with *A. thaliana* roots was studied on different media including PNM (Plant Nutrient Medium) and MS (Murashige and Skoog-Salt) by light and fluorescence microscopy without chemical treatment (Figures 1-7). With all methods, the mycelium is nicely visible in *A. thaliana* roots by optical microscope. Using different



Figure 1: *Arabidopsis thaliana* roots, and hyphae of *P. indica* under the light microscope, magnification 400x. Life microscopy. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 2: *Arabidopsis thaliana* roots, and hyphae of *P. indica* under the light microscope (A) and fluorescent channel at 450-490 nm (B), magnification 100 times. Life microscopy. 15 days after interaction of *P. indica* with *A. thaliana* on PNM medium.




Figure 3: *Arabidopsis thaliana* roots, and hyphae of *P. indica* under the light microscope (A) and fluorescent channel 450-490 nm (B), magnification 100x. Life microscopy. 15 days after co-cultivation of *P. indica* with *A. thaliana* on MS-S medium.




Figure 4: *Arabidopsis thaliana* roots, hyphae and chlamydospores of *P. indica* under the fluorescent channel of 450-490 nm (A), fluorescent and visible light (B) and light microscopy (C), magnification 100x. Life microscopy. 7 days after co-cultivation of *P. indica* with *A. thaliana* on MS-S medium.

magnifications (4x, 10x, 20x and 40x) and different wavelengths (visible and fluorescent channels) better resolution and discrimination of fungal and plant tissues could be achieved. Fluorescence microscopy is more sensitive compared to light microscopy. Additionally, different layers of the cell tissue can be scanned using fluorescence microscopy.




Figure 5: *Arabidopsis thaliana* root, and hyphae of *P. indica* under the fluorescent channel at 450-490 nm (A) and under the light microscope (B), magnification 100x. Live microscopy. 15 days after co-cultivation of *P. indica* with *A. thaliana* on MS-S medium.




Figure 6: *Arabidopsis thaliana* root, and hyphae of *P. indica* under the fluorescent channel at 450-490 nm (A) and under light microscope (B), magnification 100x. Live microscopy. 7 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 7: *Arabidopsis thaliana* root, and hyphae of *P. indica* under the fluorescence channel at 488 nm. Live microscopy. magnification 400x. 4 weeks after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.

Fluorescence microscopy

Optimal results were obtained with fuchsin acid and fluorescence microscopy using different wavelengths (Figures 8-10). The method is fast and allows visualization of *P. indica* mycelia and chlamydospores in the root environment. For detection and counting of chlamydospores, we found that staining with erythrosine and safranin using the fluorescence channel at 470 and 450-490 nm resulted in contrast between spore and roots (Figure 11A and B). With this method the mycelia are not visible and the spores can be counted without light interference. Staining with fuchsin acid/cotton blue gave comparable results. With this method, all stained parts of the roots and mycelia can be seen with visible light, but at 470 nm, only chlamydospores are visible (Figure 12).

Chemical staining

The optimal time for staining with Trypan blue (0.0001 mg/ml) is 1 minute, since longer incubations also stains root tissue. This method gave best results when only the mycelium should be visualized. It shows the interaction with root environment and allows the analysis of the mycelium networks in and around the roots. Furthermore, the strong contrast between the colored mycelia and the background is suitable for the analysis of root colonization with appropriate software (Figures 13 and 14).




Figure 8: *Arabidopsis thaliana* root, hyphae, and chlamydospores of *P. indica* under the fluorescent channel at 450-490 nm (A) and at 365 nm (B); magnification 100x. Fuchsin acid staining method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 9: *Arabidopsis thaliana* root, hyphae, and chlamydospores of *P. indica* under the fluorescent channel at 450-490 nm (A, magnification 200x) and under light microscope (B, magnification 100x). Fuchsin acid staining method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 10: *Arabidopsis thaliana* root, hyphae, and chlamydozoospores of *P. indica* under the fluorescent channel at 470 nm (A) and 450-490 nm (B), magnification 200x. Fuchsin acid staining method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 11: *Arabidopsis thaliana* root, hyphae, and chlamydozoospores of *P. indica* under the fluorescent channel at 470 nm (A), 450-490 nm (B) and under the light microscope (C), magnification 200x. Erythrosine/safranin method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 12: *Arabidopsis thaliana* root, hyphae, and chlamydozoospores of *P. indica* under the light microscope (A and C) and fluorescent channel at 470 nm (B and D), magnification 400x. Fuchsin acid / cotton blue method 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 13: *Arabidopsis thaliana* root, hyphae, and chlamydozoospores of *P. indica* light microscopy, magnification 200x. Trypan blue (0.0001 mg/ml) method (A, B, 1 min staining; C, 5 min staining). 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 14: *Arabidopsis thaliana* root, hyphae, and chlamydo spores of *P. indica* under the light microscope, magnification 200x. Trypan blue (0.01 mg/ml) method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 15: *Arabidopsis thaliana* root, hyphae, and chlamydo spores of *P. indica* under light microscope, magnification 1000x. Erythrosine/aniline blue method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.

By staining spores in the root structure we obtained optimal results with the erythrosine and aniline blue stain, since it allows visualization of both root tissue and cells (Figure 15). Other methods used were staining with janus green/methyl green/natural red, methyl green/aniline blue, naphtol blue black, fuchsin acid/trypan blue, cotton blue (Figures 14-23).

In conclusion, different staining methods allow the identification of the fungal mycelium and spores in the root environment, and partly the quantification of root colonization.

We used combinations of these staining methods to identify mutants which are impaired in the *P. indica*/*Arabidopsis* interaction and to identify plant mutants which fail to host the fungus. Besides basic research to understand the molecular mechanisms underlying root colonization and propagation of the mycelium, the comparative analysis of the staining methods may also help to analyze colonization of the roots from agriculturally important plant species.




Figure 16: *Arabidopsis thaliana* root, hyphae, and chlamydo spores of *P. indica* under the fluorescent channel at 450-490 nm (A) and under the light microscope (B), magnification 200x. Janus green/methyl green/natural red staining method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 17: *Arabidopsis thaliana* root, hyphae, and chlamydo spores of *P. indica* under the light microscope (A and C) and under the fluorescent channel at 450-490 nm (B and D), magnification 200x (A and B), 400x (C and D). Methyl green/aniline blue staining method. 15 days after interaction of *P. indica* with *A. thaliana* on PNM medium.




Figure 18: *Arabidopsis thaliana* root, hyphae, and chlamydospores of *P. indica* under the light microscope, magnification 400x. Naphtol blue black method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 19: *Arabidopsis thaliana* root, hyphae, and chlamydospores of *P. indica* under the light microscope, magnification 100x. Fuchsin acid/trypan blue method. 15 days after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.




Figure 20: *Arabidopsis thaliana* root, hyphae, and chlamydospores of *P. indica* under the light microscope, magnification 400x. Cotton blue staining method. 2 weeks after co-cultivation of *P. indica* with *A. thaliana* on PNM medium.

Acknowledgements

We like to thank Dr. Irena Sherameti, Daniela Falkenberg, Fred Ferber, Thomas Beder, Jana Rother, Franziska Eberl and Dr. Annett Eitner from Universitätsklinikum, Anatomy II.

References

- Barazani O, Benderoth M, Groten K, Kuhlemeier C, Baldwin IT (2005) *Piriformospora indica* and *Sebacina vermifera* increase growth performance at the expense of herbivore resistance in *Nicotiana attenuata*. *Oecologia*. 146:234–243.
- Camehl I, Oelmüller R. (2010) Do ethylene response factor-9 and -14 repress PR gene expression in the interaction between *Piriformospora indica* and Arabidopsis? *Plant Signal Behav*. 5:932-936.
- Camehl I, Sherameti I, Venus Y, Bethke G, Varma A, Lee J, Oelmüller R. (2010) Ethylene signalling and ethylene-targeted transcription factors are required to balance beneficial and nonbeneficial traits in the symbiosis between the endophytic fungus *Piriformospora indica* and *Arabidopsis thaliana*. *New Phytol*. 185:1062-1073.
- Deshmukh S, Hüchelhoven R, Schäfer P, Imani J, Sharma M, Waller F, Kogel KH. (2006) The root endophytic fungus *Piriformospora indica* requires host cell death for proliferation during mutualistic symbiosis with barley. *Proc Natl Acad Sci USA*. 103:18450–18457.
- Deshmukh S, Hüchelhoven R, Schäfer P, Imani J, Sharma M, Weiss M, Waller F, Kogel KH. (2006) The root endophytic fungus *Piriformospora indica* requires host cell death for proliferation during mutualistic symbiosis with barley. *Proc Natl Acad Sci USA*. 103:18450-18457.
- Glen M, Tommerup IC, Bougher NL, O'Brien PA. (2002) Are Sebacinaceae common and widespread ectomycorrhizal associates of Eucalyptus species in Australian forests? *Mycorrhiza*. 12:243–247.
- Jacobs S, Zechmann B, Molitor A, Trujillo M, Petutschnig E, Likpa V, Kogel KH, Schäfer P. (2011) Broad spectrum suppression of innate immunity is required for colonization of *Arabidopsis thaliana* roots by the fungus *Piriformospora indica*. *Plant Physiol*. 156:726-740.
- Johnson JM, Oelmüller R. (2009) Mutualism or parasitism: life in an unstable continuum. What can we learn from the mutualistic interaction between *Piriformospora indica* and *Arabidopsis thaliana*? *Endocytobiosis Cell Res*. 19:81–110.
- Kaldorf M, Koch B, Rexer KH, Kost G, Varma A. (2005) Patterns of interaction between *Populus* Esch5 and *Piriformospora indica*: A transition from mutualism to antagonism. *Plant Biol*. 7:210-218.
- Oelmüller R, Sherameti I, Tripathi S, Varma A. (2009) *Piriformospora indica*, a cultivable root endophyte with multiple biotechnological applications. *Symbiosis*. 49:1–17.
- Oelmüller R, Peškan-Berghöfer T, Shahollari B, Trebicka A, Sherameti I, Varma A. (2005) MATH domain proteins represent a novel protein family in *Arabidopsis thaliana*, and at least one member is modified in roots during the course of a plant-microbe interaction. *Physiol Plant*. 124:152–166.
- Oelmüller R, Shahollari B, Peškan-Berghöfer T, Trebicka A, Gieng PH, Sherameti I, Oudhoff M, Venus Y, Altschmied L, Varma A. (2004) Molecular analyses of the interaction between *Arabidopsis* roots and the growth-promoting fungus *Piriformospora indica*. *Endocytobiosis Cell Res*. 15:504-517.
- Peškan-Berghöfer T, Shahollari B, Gieng PH, Hehl S, Markert C, Blanke V, Varma AK, Oelmüller R. (2004) Association of *Piriformospora indica* with *Arabidopsis thaliana* roots represents a novel system to study beneficial plant-microbe interactions and involves early plant protein modifications in the endoplasmic reticulum and at the plasma membrane. *Physiol Plant*. 122:465-477.
- Pham GH, Singh A, Malla R et al. (2004) Interaction of *Piriformospora indica* with diverse microorganisms and plants. In: *Plant Surface Microbiology*, Varma A, Abbott L, Werner D, Hampp R. (eds), Springer-Verlag, Berlin, pp. 235–265.
- Rai MK, Varma A, Pandey AK. (2004) Antifungal potential of *Spiranthes calva* after inoculation of *Piriformospora indica*. *Mycoses*. 47:479-481.
- Sahay NS, Varma A. (1999) *Piriformospora indica*: a new biological hardening tool for micropropagated plants. *FEMS Microbiol Lett*. 181:297–302.
- Schäfer P, Khatabi B, Kogel KH. (2007) Root cell death and systemic

- effects of *Piriformospora indica*: a study on mutualism. *FEMS Microbiol Lett.* 275:1-7.
- Selosse MA, Bauer R, Moyersoen B. (2002a) Basal Hymenomycetes belonging to the Sebacinaceae are ectomycorrhizal on temperate deciduous trees. *New Phytologist.* 155:183-195.
- Selosse MA, Weiß M, Jany JL, Tiller A. (2002b) Communities and populations of sebacinoid basidiomycetes associated with the achlorophyllous orchid *Neottia nidus-avis* (L.) L.C.M. Rich. and neighbouring tree ectomycorrhizae. *Molecular Ecology.* 11:1183-1844.
- Selosse MA, Dubois MP, Alvarez N. (2009) Do Sebaciales commonly associate with plant roots as endophytes? *Mycol Res.* 113:1062-1069.
- Shahollari B, Vadassery J, Varma A, Oelmüller R. (2007) A leucin-rich repeat protein is required for growth promotion and enhanced seed production mediated by the endophytic fungus *Piriformospora indica* in *Arabidopsis thaliana*. *Plant J.* 50:1-13.
- Shahollari B, Varma A, Oelmüller R. (2005) Expression of a receptor kinase in *Arabidopsis* roots is stimulated by the basidiomycete *Piriformospora indica* and the protein accumulates in Triton X-100 insoluble plasma membrane microdomains. *J Plant Physiol.* 162:945-958.
- Sherameti I, Shahollari B, Venus Y, Altschmied L, Varma A, Oelmüller R. (2005) The endophytic fungus *Piriformospora indica* stimulates the expression of nitrate reductase and the starch-degrading enzyme glucan-water dikinase in tobacco and *Arabidopsis* roots through a homeodomain transcription factor which binds to a conserved motif in their promoters. *J Biol Chem.* 280:2641-2647.
- Sherameti I, Tripathi S, Varma A, Oelmüller R. (2008a) The root colonizing endophyte *Piriformospora indica* confers drought tolerance in *Arabidopsis* by stimulating the expression of drought stress-related genes in leaves. *Mol Plant-Microbe Interact.* 21:799-807.
- Sherameti I, Venus Y, Drzewiecki C, Tripathi S, Dan VM, Nitz I, Varma A, Gründler FM, Oelmüller R. (2008) PYK10, a beta-glucosidase located in the endoplasmatic reticulum, is crucial for the beneficial interaction between *Arabidopsis thaliana* and the endophytic fungus *Piriformospora indica*. *Plant J.* 54:428-439.
- Stein E, Molitor A, Kogel KH, Waller F. (2008) Systemic resistance in *Arabidopsis* conferred by the mycorrhizal fungus *Piriformospora indica* requires jasmonic acid signaling and the cytoplasmic function of NPR1. *Plant Cell Physiol.* 49:1747-1751.
- Sun C, Oelmüller R. (2010) *Piriformospora indica* hyphae and chlamydospores by Scanning Electron Microscopy. *Endocytobiosis Cell Res.* 20:34-37.
- Urban A, Weiss M, Bauer R. (2003) Ectomycorrhizas involving sebacinoid mycobionts. *Mycol Res.* 107:3-14.
- Varma A, Singh A, Sudha S, Sahay N, Sharma J, Roy A, Kumari M, Rana D, Thakran S, Deka D, Bharti K, Franken P, Hurek T, Bleichert O, Rexer KH, Kost G, Hahn A, Hock B, Maier W, Walter M, Strack D, Kranner I. (2001) *Piriformospora indica*: A cultivable mycorrhizal-like endosymbiotic fungus. In: Mycota IX. B. Bock, Springer Series, Springer-Verlag, Berlin, pp. 123-150.
- Varma A, Verma S, Sudha S, Sahay NS, Bütehorn B, Franken P. (1999) *Piriformospora indica*, a cultivable plant growth promoting root endophyte. *Appl Environ Microbiol.* 65:2741-2744.
- Verma SA, Varma A, Rexer KH, Hassel A, Kost G, Sarbhoy A, Bisen P, Bütehorn B, Franken P. (1998) *Piriformospora indica*, gen. et sp. nov., a new root-colonizing fungus. *Mycologia.* 90:898-905.
- Waller F, Achatz B, Baltruschat H, Fodor J, Becker K, Fischer M, Heier T, Hückelhoven R, Neumann C, von Wettstein D, Franken P, Kogel KH (2005) The endophytic fungus *Piriformospora indica* reprograms barley to salt-stress tolerance, disease resistance, and higher yield. *Proc Natl Acad Sci USA.* 102:13386-13391.
- Weiss M, Selosse MA, Rexer KH, Urban A, Oberwinkler F. (2004) Sebaciales: A hitherto overlooked cosm of heterobasidiomycetes with a broad mycorrhizal potential. *Mycol Res.* 108:1003-1010.
- Wells G, Huston C. (1995) High-resolution selected ion monitoring in a quadrupole ion trap mass spectrometer. *Anal Chem.* 67:3650-3655.
- Yadav V, Kumar M, Deep DK, Kumar H, Sharma R, Tripathi T, Tuteja N, Saxena AK, Johri AK. (2010) A phosphate transporter from the root endophytic fungus *Piriformospora indica* plays a role in phosphate transport to the host plant. *J Biol Chem.* 285:26532-26544.
- Zuccaro A, Basiewicz M, Zurawska M, Biedenkopf D, Kogel KH. (2009) Karyotype analysis, genome organization, and stable genetic transformation of the root colonizing fungus *Piriformospora indica*. *Fungal Genet Biol.* 46:543-550.

The Methods

Growth conditions for *Piriformospora indica* and *Arabidopsis*

Piriformospora indica was cultured as described previously (Verma et al. 1998; Peřkan-Berghöfer et al. 2004) in Petri dishes on a modified Kaefer's medium (KM). The plates were kept at room condition for 2 weeks. For co-cultivation experiments with *Arabidopsis*, *P. indica* spore solutions (1×10^6 spore/ml) were prepared from the fungus grown on KM plates. Root inoculation by *P. indica* was performed by adding 5 μ l of spore solution (1×10^6 spore/ml) to each seed on different media including MS (Murashige and Skoog) and PNM.

Life microscopy

Step 1 Plant 10 seeds of *A. thaliana* and 5 μ l of *P. indica* spores (1×10^6 s/ml) on a PNM or MS plate

Step 2 Keep at 4°C for 2 days before transfer to light (24 hour light and 23°C) for 10 days

Step 3 Observe the roots by light or fluorescent microscopy directly in the plates.

Staining with fuchsin acid

Step 1 Take samples from *Arabidopsis* root and wash in distilled water.

Step 2 Incubate the sample into fuchsin acid solution for 5 min.

Step 3 Wash the samples in distilled water for 1 min.

Step 4 Add 50 μ l GL solution on glass slide and cover with glass cover.

Step 5 Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with erythrosine/safranin

Step 1 Take samples from *Arabidopsis* root and wash in distilled water.

Step 2 Incubate the sample into erythrosine solution for 1 min.

Step 3 Stain with safranin solution for 1 min.

Step 4 Add 50 μ l GL solution on glass slide and cover with glass cover.

Step 5 Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with fuchsin acid and cotton blue

Step 1 Take samples from *Arabidopsis* root and wash in distilled water.

Step 2 Incubate the sample into fuchsin acid solution for 2 min.

Step 3 Wash the samples in distilled water for 1 min.

Step 4 Stain with cotton blue solution for 1 min.

Step 5 Washing the samples in distilled water for 1 min.

Step 6 Add 50 μ l GL solution on glass slide and cover with glass cover.

Step 7 Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with janus green/methyl green/natural red

Step 1 Take samples from *Arabidopsis* root and wash in distilled water.

Step 2 Incubate the sample into janus green and methyl green mixture solution (1:1) for 2 min.

Step 3 Wash the samples in distilled water for 1 min.

Step 4 Stain with natural red solution for 1 min.

Step 5 Add 50 μ l GL solution on glass slide and cover with glass cover.

Step 6 Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with erythrosine and aniline blue

Step 1 Take samples from *Arabidopsis* root and wash in distilled water.

Step 2 Incubate the sample into erythrosine solution for 2 min.

Step 3 Wash the samples in distilled water for 1 min.

Step 4 Stain with aniline blue solution for 1 min.

Step 5 Washing the samples in distilled water for 1 min.

Step 6 Add 50 μ l GL solution on glass slide and cover with glass cover.

Step 7 Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with trypan blue (0.0001mg/ml)

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Put the sample into trypan blue solution for 1 or 5 min.
- Step 3** Wash the samples in distilled water for 1 min.
- Step 4** Add 50 μ l GL solution on glass slide and cover with glass cover.
- Step 5** Observe the samples by light and fluorescent microscopy using different wavelengths

Staining with trypan blue (0.01 mg/ml)

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Put the sample into trypan blue solution for 1 min.
- Step 3** Wash the samples in distilled water for 1 min.
- Step 4** Add 50 μ l GL solution on glass slide and cover with glass cover.
- Step 5** Observe the samples by light and fluorescent microscopy using different wavelengths

Staining with naphthol blue black

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Incubate the sample in naphthol blue black solution for 1 min.
- Step 3** Wash the samples in distilled water for 1 min.
- Step 4** Add 50 μ l GL solution on glass slide and cover with glass cover.
- Step 5** Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with fuchsin acid and trypan blue

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Put the sample in fuchsin acid solution for 10 min.
- Step 3** Washing the samples in distilled water for 1 min.
- Step 4** Stain with trypan blue (0.0001gm/ml) solution for 3 min.
- Step 5** Washing the samples in distilled water for 1 min.
- Step 6** Adding 50 μ l GL solution on glass slide and cover with glass cover.
- Step 7** Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with methyl green and aniline blue

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Put the sample in methyl green (0.05%) solution for 10 seconds.
- Step 3** Stain with aniline blue solution for 1 min.
- Step 4** Washing the samples in distilled water for 1 min.
- Step 5** Adding 50 μ l GL solution on glass slide and cover with glass cover.
- Step 6** Observe the samples by light and fluorescent microscopy using different wavelengths.

Staining with cotton blue

- Step 1** Take samples from Arabidopsis root and wash in distilled water.
- Step 2** Put the sample in cotton blue solution for 30 seconds.
- Step 3** Put on glass slide and cover with glass cover.
- Step 4** Observe the samples by light and fluorescent microscopy using different wavelengths.

Media and Solutions

MS Medium

| | |
|-------------------------------------|---------------|
| MS-Salt (=Murashige and Skoog-Salt) | 4.4 g |
| Gel-Rite | 3 g |
| MES | 0.5 g |
| Saccharose | 13.7 g |
| H ₂ O | up to 1 liter |
| pH 5.7 – 5.8 | |

PNM Medium

| | |
|---------------------------------------|---------------|
| KNO ₃ | 5 mM |
| MgSO ₄ 1 M | 2 mM |
| Ca(NO ₃) ₂ 1 M | 2 mM |
| FeSO ₄ | 0.01 mM |
| H ₃ BO ₃ | 70 mM |
| MnCl ₂ | 14 mM |
| CuSO ₄ | 0.5 mM |
| ZnSO ₄ | 1 mM |
| Na ₂ MoO ₄ | 0.2 mM |
| CoCl ₂ | 0.01 mM |
| Gel-Rite | 3 g |
| KH ₂ PO ₄ 1 M | 2.5 ml |
| H ₂ O | up to 1 liter |
| pH 5.6 | |

Kaefer's medium

| | |
|---------------------------------|---------------|
| NaNO ₃ | 7.0 mM |
| KCl | 7.0 mM |
| MgSO ₄ | 2.1 mM |
| KH ₂ PO ₄ | 9.2 mM |
| EDTA | 0.02 mM |
| ZnSO ₄ | 0.77 mM |
| H ₃ BO ₄ | 0.18 mM |
| MnSO ₄ | 0.02 mM |
| CoCl ₂ | 0.007 mM |
| CuSO ₄ | 0.0065 mM |
| FeSO ₄ | 0.02 mM |
| ammonium molybdate | 0.001 mM |
| thiamine | 0.003 mM |
| glycine | 0.005 mM |
| nicotinic acid | 0.002 mM |
| pyridoxine | 0.0004 mM |
| glucose | 110 mM |
| peptone | 2 g/l |
| yeast extract | 1 g/l |
| casein hydrolysate | 1 g/l |
| agar | 10 g/l |
| H ₂ O | up to 1 liter |

pH 6.5

Glycerol lactic acid solution (GL)

| | |
|-----------------|-------|
| glycerol | 10 ml |
| lactic acid | 10 ml |
| distilled water | 10 ml |

Trypan blue (0.0001mg/ml) solution

| | |
|-------------|----------|
| GL solution | 1 ml |
| trypan blue | 0.0001 g |

Trypan blue 0.01mg/ml solution

| | |
|-------------|--------|
| GL solution | 1 ml |
| trypan blue | 0.01 g |

Fuchsin acid solution

| | |
|--------------|---------|
| GL solution | 1 ml |
| fuchsin acid | 0.001 g |

Erythrosine solution

| | |
|-------------|--------|
| GL solution | 1 ml |
| erythrosine | 0.01 g |

Aniline blue solution

| | |
|--------------|---------|
| GL solution | 1 ml |
| aniline blue | 0.005 g |

Naphtol blue black solution

| | |
|--------------------|----------|
| GL solution | 1 ml |
| naphtol blue black | 0.0001 g |

Safranin solution

| | |
|-------------|----------|
| LG solution | 1 ml |
| safranin | 0.0001 g |

Janus green solution

| | |
|-------------|----------|
| GL solution | 1 ml |
| janus green | 0.0001 g |

Methyl green solution

| | |
|--------------|----------|
| GL solution | 1 ml |
| methyl green | 0.0001 g |

Natural red solution

| | |
|-------------|----------|
| GL solution | 1 ml |
| natural red | 0.0001 g |

Cotton blue solution (Lactophenol cotton blue solution)

| | |
|------------------|-------|
| glycerol | 50 g |
| phenol | 25 g |
| L(+)-lactic acid | 25 g |
| cotton blue | 50 mg |